

## **68-Pin MCS Standard Connector**

Input / Output 68-Pin MCS Standard Connector

| 67 65 63 61 59 | 97531      |
|----------------|------------|
| 68 66 64 62 60 | 10 8 6 4 2 |

## MC Card 64 Channel

| Pin | 1      | GNDP (power ground)     |
|-----|--------|-------------------------|
| Pin | 2      | GNDS (signal ground)    |
| Pin | 3 - 66 | Channel 1 to 64         |
| Pin | 67     | Positive supply voltage |
| Pin | 68     | Negative supply voltage |

## MC\_Card 128 Channel Expansion

| Pin | 1      | GNDP (power ground)     |
|-----|--------|-------------------------|
| Pin | 2      | GNDS (signal ground)    |
| Pin | 3 - 66 | Channel 65 to 128       |
| Pin | 67     | Positive supply voltage |
| Pin | 68     | Negative supply voltage |

## Note:

The analog channel 64 (for a standard MEA60- or MEA120-System, or a ME-System with at least 64 channels) is available only if the digital input channel is **not** selected in the **Channel Layout** dialog box of the MC\_Rack program. Selecting the digital channel reduces the analog channels by one. For MC\_Cards with less than 128 channels, the last analog channel is not available when using the digital input channel.

Channels 61 to 64 and 125 to 128 are not used for the MEA / ME-System. Channels 61 to 63 can be used as additional analog input channels A1 to A3 via the additional BNC connectors. ME-Systems may not support all 64 channels, depending on the version of the MC\_Card. For example, a ME16-System supports only channels 1 to 16.

Operating an USB supported system without the data acquisition board MC\_Card, the digital channel(s) are additional available without reducing the analog channels.